

## Precision Balance



### **SPECIFICATIONS:**

Model	PG-3000	PG-5000
Make	Galaxy	Galaxy
Capacity	3000 g	5000 g
Resolution	0.01 g	0.01 g
Calibration	External	External
Pan Size	Ø 130 mm	Ø 130 mm
Display	LCD	LCD
Response Time	3-4 secs	3-4 secs
Units	g, ct, oz	g, ct, oz
Battery Backup	In-built	In-built
Dimensions	205 x 272 x 163 mm	
Weight	3 kg	3 kg

**ACCESSORIES:** • Standard Weight • Adapter