

## Analytical Balance

1mg (EMFC) Professional Balance



### FEATURES HIGHLIGHTS:

- Electro Magnetic Force Compensation
- Under-hook Weighing
- Internal Auto Calibration
- LCD Display
- Multiple Weighing Units
- Tare
- RS232C Interface
- Counting Weighing Model
- Density Determination (Manual Calculation)

## SPECIFICATIONS:

Model	MA-201	MA-501	MA-1001
Capacity	220g	520g	1000g
Make	Galaxy		
Readability	0.001g (1mg)		
Linearity	± 0.002g		
Minimum Load	0.1 g		
Display	LCD		
Pan Size	Ø 90 mm	Ø 90 mm	Ø 120 mm
Response Time	4 ~ 5 sec		
Calibration	Internal		
Tare Range	Full		
Operating Temperature	15 to 35°C		
Operating Humidity	40~ 80% R.H		
Warm Up time	30 min		
Units	mg, g, ct, ozt, oz, GN, dwt, lb		
Key Sound	Yes		
Weighing Chamber Height	195 mm		
Baud Rate	600, 1200, 2400, 4800, 9600		
Class	I		
Interface	RS232C		
Date Output	Command and continuous output		
Power Supply	AC 220V / 50 Hz		
Dimensions & Weight	220 x 360 x 345 mm & 10 kg		

### STANDARD

- Draft Ring
- Power Cord

### ACCESSORIES:

- Pan
- Instruction Manual