


Precision Balance

300g x 1mg (0.001g)



MODEL: PG-300, PG-500

MODEL: PG-300S

SPECIFICATIONS:

Model	PG-300	PG-500	PG-300S
Make	Galaxy	Galaxy	Galaxy
Capacity	300 g	500 g	300 g
Readability	0.001g (1mg)	0.001g (1mg)	0.001g (1mg)
Repeatability	0.001 g	0.001 g	0.001 g
Linearity	0.002 g	0.002 g	0.002 g
Pan Size	Ø 90 mm	Ø 90 mm	Ø 90 mm
Response Time	3 - 4 sec	3 - 4 sec	3 - 4 sec
Calibration	External	External	External
Display	Backlight LCD	Backlight LCD	LED
Tare Range	Full	Full	Full
Operating Temperature	10 to 40°C	10 to 40°C	10 to 40°C
Power Supply	AC 220V / 50Hz	AC 220V / 50Hz	AC 220V / 50Hz
Dimensions	325 x 300 x 225	325 x 300 x 225	190 x 240 x 100
Weight	4 kg	4 kg	3 kg

ACCESSORIES: • Standard Weight • Adapter