

Precision Balance

(Economy)



This is a stable & accurate 0.01g (10mg) precision weighing balance. With inbuilt battery & standard weights it is an ideal choice for every workplace.

ACCESSORIES: • Standard Weight • Adapter

SPECIFICATIONS:

Model	PG-600	PG-1000
Make	Galaxy	
Capacity	600 g	1000 g
Resolution	0.01 g	
Calibration	External	
Pan Size	Ø 128 mm	
Display	LED	
Response Time	3-4 secs	
Units	g, ct, oz	
Battery Backup	In-built	
Dimensions	190 x 240 x 100	
Weight	3 kg	