


Precision Balance



SALIENT FEATURES:

- Green LED Display
- Front & Back display
- With Rs232 data output
- Removable glass windshield
- Low battery indicator
- In-built rechargeable battery for continuous use
- Piece counting & percentage function
- Multiple weighing units
- Auto zero tracking

SPECIFICATIONS:

| Model | PG-620 | PG-1050 |
|-----------------|----------------------------------|---------|
| Make | Galaxy | |
| Capacity | 620 g | 1050 g |
| Resolution | 0.01 g | |
| Repeatability | ± 0.01 | |
| Linearity | ± 0.02 | |
| Calibration | External Calibration | |
| Pan Size | Ø 120 mm | |
| Tare Range | Full | |
| Operating Temp. | 10 to 40°C | |
| Display | Green LED (Front & back display) | |
| Response Time | 3-4 secs | |
| Units | g, ct, oz, dwt | |
| Battery Backup | In-built | |
| Dimensions | 190 x 240 x 100 | |
| Weight | 3 kg | |

ACCESSORIES:

- Windshield Set
- Pan
- Instruction Manual
- Adapter