

pH Meter Economy

MODEL: PH-3M

Microprocessor based



SALIENT FEATURES:

- Manual Temperature Compensation
- Large LCD display with blue backlight.
- Simultaneously displays pH & Temperature.
- Smart functions for calibration, data storage
- Automatically displays the percentage of slope after calibration is over
- Equipped with electrode holder, calibration powder & electrode

SPECIFICATIONS:

Model	PH-3M
Make	Galaxy
Range	0 to 14.00 pH
pH Resolution	0.1 pH/ 0.01 pH
Accuracy	± 0.01 pH
Calibration Points	up to 3 points
pH Buffer options	USA (pH 4.01, 7.00, 10.01)
mV Range	± 1999 mV
Accuracy	± 1% FS ±1 digit
Resolution	1 mV
Temperature Range	0 to 100°C
Accuracy	± 0.5°C
Resolution	0.1°C
Temperature Compensation	0 to 100°C, Manual
Connector	BNC
Display	Backlight LCD
Power	DC9v, AC220v, 50Hz
Dimensions & Weight	160 x 190 x 170 mm, 750g

STANDARD

ACCESSORIES:

- pH Electrode
- pH Buffer
- Electrode Stand
- Adapter
- Instruction Manual