

Accessories for Moisture Analyzer

Pans for Moisture Analyzer

MODEL: PMB-60



These aluminum disposable pans are suitable for all moisture balances.

It's 90 mm diameter and 8 mm depth serves well for all machines.

It's oil free surface results in accurate moisture analysis.

80 pcs/box

Glass Fiber Pads for Moisture Analyzer

MODEL: GFP-60



These glass fiber pads are suitable for all moisture balances.

Dimension: 90 mm (ø)

They increase the evaporation surface of a sample and reduce measurement time. The filters are particularly useful when testing liquids, semi-liquids and substances in case of which crust formation takes place.

60pcs/box

Test Standard for Moisture Analyzer

MODEL: TSMA-1



The moisture content rate of sodium tartrate dihydrate is 15.66%. Since its moisture content is theoretically known, it can be used as a test sample to confirm the performance of analyzers.

Sodium Tartarate Dihydrate - $C_4H_8Na_2O_8$

1 Test - 20g/pack