

## Accessories for Electrochemical Testers

### pH Electrode



Model	PE-4
Range	0~14 pH
Temp. Range	0~60°C
Reference	AgCl
Body & Connector	Plastic, BNC
Dimensions	12(∅) x 150(L)mm

### pH Electrode (suitable for 5 point)



Model	PE-5
Range	0~14 pH
Temp. Range	0~100°C
Reference	AgCl
Body & Connector	Glass, BNC
Dimensions	12(∅) x 150(L)mm

### TDS/Conductivity Cell



Model	CON-1	CON-10	CON-0.1	CON-1-3P
Range	10 μS/cm to 20 mS/cm	0.1 to 200 mS/cm	0 to 100 μS/cm	10 μS/cm to 20 mS/cm
Cell Constant	K=1	K=10	K=0.1	K=1
Material	Platinum sensor, glass body			
Temp. Range	0 to 80°C	0 to 80°C	0 to 80°C	0 to 80°C
Connector	6-pin mini-DIN	6-pin mini-DIN	6-pin mini-DIN	3-pin mini-DIN
Dimensions	12(∅) x 150(L)mm	12(∅) x 150(L)mm	12(∅) x 150(L)mm	12(∅) x 150(L)mm

### Dissolved Oxygen Probe



Model	DO-1
Range	0~20mg/L
Sensor	Polarographic
Temp.	0~50°C
Connector	6-pin mini-DIN
Dimensions	12(∅) x 150(L)mm

### Temperature Probe



Model	TP-10K
Temp. Range	0 to 105°C
Sensor	NTC 10 KΩ
Body	S. steel and plastic
Connector	3.5 mm jack plug
Dimensions	120 × 3(∅) mm