

## Precision Balance



**PG-10000**



**PG-3000, PG-5000**

### SPECIFICATIONS:

Model	PG-10000	PG-3000	PG-5000
Make	Galaxy	Galaxy	Galaxy
Capacity	10000 g	3000 g	5000 g
Resolution	0.1 g	0.01 g	0.01 g
Repeatability	± 0.1 g	± 0.01 g	± 0.01 g
Linearity	± 0.2 g	± 0.02 g	± 0.02 g
Calibration	External	External	External
Pan Size	145 x 175 mm	Ø 130 mm	Ø 130 mm
Display	LCD	LCD	LCD
Response Time	3-4 secs	3-4 secs	3-4 secs
Units	g, ct, oz	g, ct, oz	g, ct, oz
Battery Backup	In-built	In-built	In-built
Dimensions	205 x 272 x 163 mm		
Weight	3 kg	3 kg	3 kg

**ACCESSORIES:** • Standard Weight • Adapter