

Digital Portable Refractometer

MODEL: PR-90



SALIENT FEATURES:

- Automatic 1 touch measurement
 - Can be calibrated to 0 using pure water
 - Automatic temperature compensation
 - 2 in 1: Brix & nD Combo
- Can change scale by single touch

SPECIFICATIONS:

Model	PR-90
Make	Galaxy
Measurement range	Brix: 0~90 % nD: 1.3330 ~ 1.5177
Resolution	Brix: 0.1 %, nD: 0.0001
Accuracy	Brix: $\pm 0.1\%$, nD: ± 0.0003
Measurement range (Temp.)	0 ~ 80°C
Resolution	0.1°C
Accuracy	$\pm 1^\circ\text{C}$
Light Source	Yellow LED
Measurement Time	1.5 sec
Sample Volume	300 μL
Power Source	2 x 1.5v(AAA)
Dimensions & Weight	130 x 55 x 25 mm - 0.2 kg

STANDARD

ACCESSORIES:

- Main Unit
- Pasteur Pipette
- Screwdriver
- Batteries
- Instruction Manual