


## Low Viscosity Adapter

for Viscometers

**MODEL: LVA-30**



SALIENT FEATURES:

- Ideal for low viscosity materials
- 20 to 30 mL sample is needed for test
- Spindle No. included with adapter
- Especially made for materials less than 20 mpa.s viscosity
- Viscosity Measurement Range: 1 to 2000 mpa.s

STANDARD ACCESSORIES:

- Joint Road
- Hook Rod
- Frame
- Spindle No. 0
- Pipe

## Circulating Jacketed Container

for Viscometers

**MODEL: CJC-200**



- Useful for maintaining samples at a constant temperature
- Connect this to a circulating water/oil bath and maintain your sample at constant temperature
- Manufactured from heat resistant borosilicate glass
- Flat surface and free of bubbles

## Temperature Sensor for Viscometers

**TMP-2**

compatible for GV-60

**TMP-3**

compatible for GV-100 & GV-200

