



Lab grade series is ideal for applications where operation time is shorter than 30 minutes. These ultrasonic cleaners can be used for a maximum of 30 minutes at a stretch, after 30 minutes of operation, 15 minutes of rest time is required before it can be used again.

Lab Grade (for short operating time)

Model	UC-1002	UC-1003	UC-1004	UC-1006	UC-1010	UC-1015	UC-1022	UC-1030
Theoretical Capacity	2.0L	3.2L	4.5L	6.5L	10.0L	15.0L	22.0L	30.0L
Actual water filling	1.6 L	2.5 L	3.5 L	5 L	8 L	12 L	18 L	27 L
Tank Size (mm)	150*135*100	240*135*100	300*150*100	300*150*150	300*240*150	330*300*150	500*300*150	500*300*200
Overall Size (mm)	180*165*205	270*185*230	330*180*240	375*200*275	330*270*310	360*330*310	530*330*310	530*325*325
Ultrasonic Power	60W	120W	180W	180W	240W	360W	480W	600W
Heating Power	100W	100W	200W	200W	200W	300W	500W	500W
Tank Drainage Valve	No	No	No	Yes	Yes	Yes	Yes	Yes
Time Setting	1-30 mins	1-30 mins	1-30 mins	1-30 mins	1-30 mins	1-30 mins	1-30 mins	1-30 mins
Heating Temp.	0-60°C	0-60°C	0-60°C	0-60°C	0-60°C	0-60°C	0-60°C	0-60°C
Display	LED	LED	LED	LED	LED	LED	LED	LED
Frequency	40kHz	40kHz	40kHz	40kHz	40kHz	40kHz	40kHz	40kHz
Lid Material	SS 204	SS 204	SS 204	SS 204	SS 204	SS 204	SS 204	SS 204
Inner Tank Material	SS 304	SS 304	SS 304	SS 304	SS 304	SS 304	SS 304	SS 304
Power Supply	AC 220V	AC 220V	AC 220V	AC 220V	AC 220V	AC 220V	AC 220V	AC 220V
Weight	2.4kg	3.6kg	4.9kg	5.4kg	7.5kg	10kg	11kg	18kg

STANDARD

- Main Instrument
- Basket
- Instruction Manual

ACCESSORIES:

- Lid
- Power Cord