

Bacti-Cinerator Loop Sterilizer

MODEL: BCS-150



SALIENT FEATURES:

- Compact size and max sterility without splatter across the work surface area
- Sterilizes needles, loops, and culture tube mouths
- Flash sterilization in 5 to 7 seconds at optimum sterilizing temperature
- Ceramic funnel tube enclosed in stainless-steel perforated guard and stainless-steel support stand

SPECIFICATIONS:

Model	BCS-150
Make	Galaxy
Central highest temperature	825°C ±50°C
Max Diameter of object to be sterilized	Ø 14mm
Length of heating zone	150 mm
Modes	Low Temperature Mode 400°C
	High Temperature Mode 825°C
Heating Time 400°C	10 min
Heating Time 825°C	20 min
Sterilization Time	5 to 7 seconds
Power Supply	AC 220V / 50 Hz
Dimensions & Weight	96 × 192 × 195 mm & 1.3 kg

STANDARD

ACCESSORIES:

- Power Cord
- Instruction Manual