

## Roller Mixer

**MODEL: GRM-6**



### SALIENT FEATURES:

- Rocking and rolling motion for complete mixing
- Digital control of time and speed
- Six convenient rollers
- Roller length 28 cm
- Suitable for various tube diameters
- Stable functioning and very low noise

### SPECIFICATIONS:

Model	GRM-6
Make	Galaxy
Motion	Rocking and Rolling
Speed Range	20 to 80 RPM
Rollers	6 rollers
Max Capacity	10 * 50ml Centrifuge tube / 0.5kg
Timing Range	1min to 23h 59mins
Compatible Tube Dia	Ø6 to 34mm
Roller Length	28 cm
Power	12W
Power Supply	AC 220V / 50 Hz
Dimensions & Weight	445 × 260 × 103 mm & 3 kg

### STANDARD

### ACCESSORIES:

- Adapter
- Instruction Manual