

Tube Rotator Mixer



MODEL: GTR-H



MODEL: GTR-D

SALIENT FEATURES:

- Compatible with a variety of tube clamps, suitable for 1.5 mL to 50 mL tubes
- Digital control of time and speed
- Removable tube holder enables easy cleaning and maintenance in case of sample spillage
- Adjustable rotor tilt angle from 0° to 90°
- Gentle mixing action minimizes the risk of cell rupture, ideal for delicate biological samples
- Stable functioning and very low noise
- Designed for continuous, stable mixing at a constant speed over long durations.

SPECIFICATIONS:

Model	GTR-H	GTR-D
Make	Galaxy	Galaxy
Motion	360° Rotating	360° Rotating
Speed Range	20 to 80 RPM	20 to 80 RPM
Tilt Angle	0 to 90°	0 to 90°
Max Capacity	0.8 kg	0.8 kg
Timing Range	1min to 99h 59mins	1min to 23h 59mins
Compatible Tube Dia	1.5ml to 50ml	1.5ml to 50ml
Power	12W	12W
Power Supply	AC 220V / 50 Hz	AC 220V / 50 Hz
Dimensions & Weight	415 × 110 × 140 mm & 3.5 kg	230 × 250 × 260 mm & 2 kg

STANDARD

- Adapter

ACCESSORIES:

- Instruction Manual