

Ice Flaker Machine

IFM-300



SALIENT FEATURES:

- Fully Automatic Microprocessor control.
- Stainless Steel construction.
- Continuous ice flakes output
- CFC free compressor
- Low water level indication
- Overload protection
- Noiseless operation
- Compact Design and low maintenance

SPECIFICATIONS:

Model	IFM-300
Make	Galaxy
Capacity of Ice Flaker	300 kg / 24 hr
Storage Capacity	200 kg
Storage Bin	ABS Plastic
Cooling Method	Air (inside R404A refrigerant)
Ambient Temperature	Below 25°C for optimum working
Construction	Stainless Steel
Noise Level	≤52 dB
Water Supply	3/4" BSP
Refrigerant	R404A
Power Supply	AC 220V / 50 Hz
Input Power	1360 W
Dimensions (mm)	560 x 800 x 1700 mm
Gross Weight	105 kgs