

Analytical Balance

220g x 0.1mg

MODEL: MA-221



FEATURES HIGHLIGHTS:

- Electro Magnetic Force Compensation
- Under-hook Weighing
- Internal Auto Calibration
- LCD Display
- Multiple Weighing Units
- Tare
- RS232C Interface
- Counting Weighing Model
- Density Determination (Manual calculation)

SPECIFICATIONS:

Model	MA-221
Capacity	220g
Make	Galaxy
Readability	0.0001g (0.1mg)
Linearity	± 0.0003g
Minimum Load	0.01g
Display	LCD
Pan Size	Ø 90 mm
Response Time	4-6 sec
Calibration	Internal Calibration
Tare Range	Full
Operating Temperature	15 to 30°C
Operating Humidity	50-70% R.H
Warm Up time	30 min
Units	g, ct, ozt, oz, GN, dwt, lb, mg
Key Sound	Yes
Weighing Chamber Height	210 mm
Baud Rate	1200, 2400, 4800, 9600
Class	I
Interface	RS-232C
Date Output	Command & Continuous
Power Supply	AC 220V / 50 Hz
Dimensions & Weight	220 x 360 x 345 mm & 9.5 kg

STANDARD

- Draft Ring
- Instruction Manual

ACCESSORIES:

- Pan
- Power Cord