

Moisture Analyser

MODEL: MB-6900



SALIENT FEATURES:

- Intuitive touch screen
- Halogen Light heating
- Temperature & time can be set
- Stores 20 data records
- Automatically displays the percentage of slop
- Displays moisture %, dry %, dry weight
- 5 inch touchscreen
- 3 ending modes and 3 drying modes

SPECIFICATIONS:

Model	MB-6900
Make	Galaxy
Moisture Range	0~100%
Moisture Resolution & Accuracy	R:0.001%, A: ±0.05%
Capacity	110g
Readability	0.001 g (1 mg)
Repeatability & Linearity	± 0.001g & ± 0.002g
Minimum Load	0.01g
Display & Pan Size	5" Touchscreen, Ø 100 mm
Sensor & Calibration	Loadcell, External
Heating Source	Halogen Lamp
Drying Modes	Standard, Fast & Soft
Ending Modes	Manual, Timed, Automatic
Temp. Setting Range	60 to 200°C
Time Setting	1 to 99 mins
No. of Storage	20 data sets
Interface	Rs232
Power Supply	AC 220V / 50 Hz
Dimensions & Weight	350 × 215 × 190 mm & 4.5 kg

STANDARD

ACCESSORIES:

- 20 x Pan
- Pan Holder
- Pan Lifter
- Adapter
- Instruction Manual