

## Mini Centrifuge

Fixed speed - 6000 RPM

**MODEL: MCX-6F**



### SALIENT FEATURES:

- Brushless motor for long life.
- Can use 2/1.5ml micro centrifuge tube as well as PCR tubes.
- Automatic motor shut down on opening the cover.
- Compact and easy-to-use design.

### SPECIFICATIONS:

| Model                      | MCX-6F                        |
|----------------------------|-------------------------------|
| Make                       | Galaxy                        |
| Speed                      | Fixed 6000 RPM                |
| Relative centrifugal force | 3300 g                        |
| Timer                      | 1 to 99 minutes               |
| Motor Type                 | Brushless DC Motor            |
| Motor Power                | 48 W                          |
| Noise                      | 55 dBA                        |
| Rotor                      | 1 All-in-one rotor            |
| Tubes which can be used    | MCT 6*2/1.5: mL, PCR 16*0.2mL |
| Power Supply               | DC 24v/2A, AC 220v / 50Hz     |
| Dimensions & Weight        | 153 x 155 x 112 mm & 1 kg     |

### STANDARD

### ACCESSORIES:

- Standard Rotor
- Power Adapter
- Instruction Manual