

## Vortex Mixer

**MODEL: VM-303**



### SALIENT FEATURES:

- Suitable for 0.2ml to 50ml micro tubes  
Also suitable for containers less than 108mm diameter.
- Continuous & touch operation
- Brush-less DC motor for maintenance free long life
- Built in counter balance for steady & safe operation

### SPECIFICATIONS:

Model	VM-303
Make	Galaxy
Rotating Speed	300 to 4200 RPM
Shaking Movement	Orbital
Orbital Diameter	4 mm
Motor Power	48 W
Run Type	Continuous / Touch / Speed
Pulse Mode	Adjustable
Max. Load Capacity	500 g
Ambient Temp. Range	5 to 40°C
Ambient Humidity	80%
Power Supply	AC 220V / 50 Hz
Dimensions & Weight	120 × 170 × 160 mm & 2.2 kg

### STANDARD

### ACCESSORIES:

- Standard Head - A
- Flat Head - B
- Foam Module - E
- Adapter
- Instruction Manual

