


Standard Weight Set as per OIML R111




Configuration: 200g x 2, 100g x 1, 50g x 1, 20g x 2,
10g x 1, 5g x 1, 2g x 2, 1g x 1,
500mg x 1, 200mg x 2, 100mg x 1,
50mg x 1, 20mg x 2, 10mg x 1,
5mg x 1, 2mg x 2, 1mg x 1

SHAPES OF MILLIGRAM WEIGHTS:



1 mg
10 mg
100 mg

Triangle
1 segment



2 mg
20 mg
200 mg

Square
2 segments



5 mg
50 mg
500 mg

Pentagon
5 segments

SPECIFICATIONS:

Model	WS-E1-200	WS-E2-200	WS-F1-200
Class	E1	E2	F1
Standard	OIML R111		
Material	Austenitic Stainless Steel	316 Stainless Steel	304 Stainless Steel
Construction	Solid single piece fabrication		1 or 2 pieces
Cavity	No cavity adjustment		May contain
Density	8.00g/cm ³	7.96g/cm ³	7.9g/cm ³
Magnetic Susceptibility	<0.005	<0.005	<0.05
Intensity of polarization	<2.5 μ T	<8 μ T	<25 μ T
Shape	Cylindrical with knob		
Finish	Mirror Polish		
Packing	ABS & Aluminum box		
Accessories	Tweezer, Cloth, Brush & Gloves		Tweezer & Cloth