

Small Sample Adapter

for Viscometers

Model: SSA-LV

- Spindle No. 18, 25, 31, 34
- Suitable for LV Series Viscometer
- Also suitable for use with High Temperature Cell
- Require just 7 to 11 ml of sample

Model: SSA-RV

- Spindle No. 21, 27, 28, 29
- Suitable for RV/HV Series Viscometer
- Also suitable for use with High Temperature Cell
- Require just 8 to 13 ml of sample



Circulating Jacketed Container

for Viscometers

MODEL: CJC-200



- Useful for maintaining samples at a constant temperature
- Connect this to a circulating water/oil bath and maintain your sample at constant temperature
- Manufactured from heat resistant borosilicate glass
- Flat surface and free of bubbles

High Temperature Cell

for Viscometers

MODEL: VTC-1



- Useful for heating up and maintaining samples at a constant temperature
- PID controller ensures accurate temperature control

Temperature Sensor for Viscometers

TMP-2

compatible for GV-60



TMP-3

compatible for GV-100 & GV-200

