


Standard Weight Set as per OIML R111



Configuration: 25 pieces


2000g x 2, 1000g x 1, 500g x 1, 200g x 2,
100g x 1, 50g x 1, 20g x 2, 10g x 1, 5g x 1,
2g x 2, 1g x 1, 500mg x 1, 200mg x 2,
100mg x 1, 50mg x 1, 20mg x 2, 10mg x 1,
5mg x 1, 2mg x 2, 1mg x 1

SHAPES OF MILIGRAM WEIGHTS:



1 mg
10 mg
100 mg

Triangle
1 segment



2 mg
20 mg
200 mg

Square
2 segments



5 mg
50 mg
500 mg

Pentagon
5 segments

SPECIFICATIONS:

Model	WS-F1-2000
Class	F1
Standard	OIML R111
Material	304 Stainless Steel
Construction	1 or 2 pieces
Cavity	May contain
Density	7.9g/cm ³
Magnetic Susceptibility	<0.05
Intensity of polarization	<25 μT
Shape	Cylinderical with knob
Finish	Mirror Polish
Packing	ABS & Aluminum box
Accessories	Tweezer & Cloth