


Precision Balance



Salient Features

- In-built battery backup
- Backlight display for easy viewing
- Compact and sturdy
- Stable reading

SPECIFICATIONS:

Model	PG-300M	PG-500M
Make	Acutek	Acutek
Capacity	300 g	500 g
Readability	0.001g (1mg)	0.001g (1mg)
Repeatability	0.001 g	0.001 g
Linearity	0.002 g	0.002 g
Pan Size	Ø 90 mm	Ø 90 mm
Response Time	3 - 4 sec	3 - 4 sec
Calibration	External	External
Display	Backlight LCD	Backlight LCD
Tare Range	Full	Full
Operating Temperature	10 to 40°C	10 to 40°C
Power Supply	AC 220V / 50Hz	AC 220V / 50Hz
Dimensions	325 x 300 x 225	325 x 300 x 225
Weight	4 kg	4 kg

ACCESSORIES: • Standard Weight • Adapter